



Mobility

3 & 4 Wheeled Walkers	2-4
Wheelchairs	5
Walking Sticks	6-9
Walking Frames	10
Walker Accessories	11-13
Handy Bar	14
Ramps	15

1 Tri Wheels

High quality 3 wheeled walker made from strong gauge steel. Folds for easy storage and transportation. Powder coated frame. Orthopaedic handles with hook, lock down brakes. Optional baskets and tray. Height adjustable handles. Available in three frame colours. Solid tyres.

Lowest Handle Height	863mm (34")
Highest Handle Height	940mm (37")
Width at Widest Point of Tri wheel	610mm (24")
Carrying Weight	6.5kg (14lbs)

PR30270/30-C Chrome
 PR30273/30-BL Blue
 PR30273/30-BU Burgundy



2 Aluminium Tri Wheel Walker

High quality 3 wheeled walker made from aluminium. Folds for easy storage and transportation. Powder coated frame. Orthopaedic handles with hook lock down brakes. Optional basket and tray. Height adjustable handles. Available in Blue, Burgundy and Silver. With solid tyres.

PR30274/30-AL Blue
 PR30274/30-AL/SL Silver
 PR30274/30-AL/BU Burgundy

Lowest Handle Height	736mm (29")
Highest Handle Height	940mm (37")
Width at Widest Point of Tri wheel	585mm (23")
Carrying Weight	4.2kg (9lbs)



3 Folding 3 Wheel Walker - with push down brakes

High quality 3 wheeled walker. Folds for easy storage and transportation. Orthopaedic round handles with push down brakes. Height adjustable handles. Solid tyres. Chrome.

Lowest Handle Height	863mm (34")
Highest Handle Height	940mm (37")
Width at Widest Point of Tri wheel	630mm (25")
Carrying Weight	6.5kg (14lbs)

PR30270/40



1

1 Aluminium Rollator - Blue

A light weight aluminium rollator that is easily manoeuvrable, prompting independence inside and outside of the home. Looped brakes can be locked to give the client additional confidence when resting in a seated or standing position. Includes a cane holder attachment. Folds easily for storage and transportation. Height adjustable handles. Available in a metallic blue frame. Includes free cane holder.

Seat Height	584mm (23")
Handle Height	841-951mm (33"-37½")
Width Between Handles	458mm (18")
External Width	635mm (25")
Weight	7.6kg (16.7 lbs)
Wheel Size	204mm (8")

PR30283/30-BL
PR30283/30-2 Replacement Basket

Max user weight  



2

2 Aluminium Rollator - Burgundy

A lightweight aluminium rollator, that is easily manoeuvrable, promoting independence inside and outside the home. Hook brakes. Includes seat with comfortable shaped back rest and removable basket. Folds easily for storage and transportation. Colour Metallic Burgundy.

Seat Height	584mm (23")
Handle Height	830-940mm (32½"-37")
Width Between Handles	458mm (18")
External Width	635mm (25")
Weight	6.5kg (14 lbs)
Wheel Size	165mm (6½")

PR30383/36-BU

Max user weight  



3

3 Aluminium Rollator - Champagne

This stylish champagne coloured rollator has a height adjustable frame that allows the seat to be height adjustable to suit the user. The handles and frame adjust independently, giving an excellent range of height options, supplied with four large castors for easy manoeuvrability and cable brakes giving superior control, especially when used outdoors. The frame folds easily for transport and storage.

Seat Height	570-620mm (22½"-24½")
Handle Height	825-965mm (32½"-38")
Width Between Handles	440mm (17½")
External Width	540mm (21½")
Weight	6.6kg 6.5kg (14½ lbs)
Wheel Size	178mm (7")

PR30279

Max user weight  

1 Lets Go Indoor Rollator

The Lets Go Indoor Rollator is an entirely new and practical walking aid. It is light, flexible and beautifully designed for indoor use in the home, hospital or care home.

It is comfortable to use and easy to manoeuvre over carpet and household thresholds, and it is perfect for getting into narrow spaces such as toilets and bathrooms.

The hand operated brakes are easy to apply and can be used as parking brakes. The handle is soft and comfortable and can be adjusted to suit the height of the user. The brakes are also linked together so that both brakes can be operated with just one hand.

The Lets Go comes with a useful carry bag and clip on tray, which is ideal for transporting drinks, plates or any other small light weight items from room to room.

The Lets Go can be neatly folded away when not in use or for transporting in a car.

A useful accessory in the hospital environment is the drip holder. This optional extra is easily attached to the frame and allows the user to walk with a drip and a walking aid comfortably.

Aluminium finish with stylish beech appearance on the rear legs.

Handle Height	790-900mm (25½"-31")
Tray Height	680mm (27")
Length	630mm (25")
Width	520mm (20½")
Width between Wheel at Front	400mm (16")
Width between Wheel at Back	480mm (19")
Wheel Diameter	150mm (6")
Carrying Weight	6kg (13lbs)

- PR30285 Let Go Indoor Rollator
- PR30285/T Replacement Tray
- PR30285/FB Replacement Fabric Bag
- PR30285/DB Drip Bag Holder

Max user weight  



2 Deluxe Rollator

This Deluxe Rollator is constructed with anodised lightweight aluminium. Its novel folding mechanism along with patented quick release front wheels ensure that when not in use it can be stored and transported in the smallest of areas. It has "one touch" locking loop braking system, and four 8" PU tyres, with the front two as swivel wheels.

A mesh cloth back rest and seat provide comfort when sitting, and a removable shopping bag is also included.

Colour: Metallic Blue

Handle Height	790-915mm (26"-31")
Seat Height	533-538mm (21"-21")
Length	660-670mm (26"-26")
Width between Wheel at Front	652-662mm (25½"-26")
Wheel Diameter	200mm (8")
Carrying Weight	7.94kg (17.5 lbs)

PR30280

Max user weight  





Folding Adjustable Walking Sticks

These folding adjustable walking sticks with “T” shaped handle are lightweight and easy to use, they fold away neatly for easy storage when not in use. All canes have a maximum user weight of 114 kg (18 stone). Ferrule size for these canes is 19mm (¾”)

Code	Colour	Size
PR30121/I	Birds Eye Maple	838-940mm (33-37")
PR30121/J	Natural	838-940mm (33-37")
PR30121-F/ICE	Blue Ice	838-940mm (33-37")
PR30121-F/ICE/S	Blue Ice	737-838mm (29-33")
PR30120-F/FRL	Floral	787-889mm (31-35")
PR30120-F/FRL/S	Floral	737-838mm (29-33")
PR30120-BK	Black	787-889mm (31-35")
PR30120-BK/S	Black	737-838mm (29-33")
PR30120-BZ	Bronze	787-889mm (31-35")
PR30121-BK	Black	838-940mm (33-37")
PR30121-BZ	Bronze	838-940mm (33-37")
PR30120/BF	Burgundy Flower	838-940mm (33-37")
PR30120/BF/S	Burgundy Flower	737-838mm (29-33")
PR30120/BP	Blue Petal	838-940mm (33-37")
PR30120/BP/S	Blue Petal	737-838mm (29-33")
PR30120/T	Tartan	838-940mm (33-37")
PR30120/T/S	Tartan	737-838mm (29-33")



Max user weight   All sticks & canes on this page



1 Arthritis Grip Cane Adjustable - Black

Arthritis Grip Cane for users that need full palm support. Especially suitable for Rheumatoid Arthritis sufferers. Adjustable in height 762-991mm (30-39"). Ferrule size for these canes is 19mm (¾").


PR30107-BK Right Handed
PR30108-BK Left Handed



2 Arthritis Grip Cane - Adjustable - Folding - Brown

Arthritis grip cane for users that need full palm support. Especially suitable for Rheumatoid Arthritis sufferers. Adjustable in height 838-940mm (33-37"). Cane folds for easy storage. Ferrule size for these canes is 19mm (¾").

PR30107F-BR Right Handed
PR30108F-BR Left Handed



3 Adjustable Quad Cane

Height adjustable quad cane with curved neck, which places the users weight directly over the shaft for safety. The handle is shaped to fit the centre of the palm and the front tapers for the fingers. Steel construction with graphite coloured shaft and black powder coating finish at the base for durability. Extra grip rubber ferrules on all four legs. Height adjusts from 750-930mm (29.5-36.5"). Base dimension 230 x 240mm (9 x 9.5").

PR30148/SB



4 Comfort Grip Cane - Adjustable - Folding - Brown

Comfort grip cane for users that need gentle palm support. Adjustable in height 838-940mm (33-37"). Cane folds for easy storage. Ferrule size for these canes is 19mm (¾").


PR30107FS-BR Right Handed
PR30108FS-BR Left Handed



5 Comfort Grip Cane Adjustable - Small Handle - Black

Comfort Grip Cane for users that need gentle palm support. Adjustable in height 762-991mm (30-39"). Ferrule size for these canes is 19mm (¾").

PR30107S-BK Right Handed
PR30108S-BK Left Handed



6 Adjustable Aluminium Cane

Adjustable aluminium cane with plastic York handle. Height adjustable from 787-991mm (31-39"). Bronze in colour. Ferrule size for this canes is 22mm (7/8").

PR30130-BZ

7 Non-Folding Aluminium Adjustable Cane - Wood Grain Mahogany

T shape wooden handle. Light and easy to use, height adjusts from 787-991mm (31-39"). Mahogany. Ferrule size for this canes is 22mm (7/8").

PR30090-E

Max user weight All sticks & canes on this page

1 Wooden Walking Canes

These wooden walking canes are beautifully crafted from beech, and finished with a Derby hook shaped handle and rubber ferrule. They all measure 36" (915mm) in length. Ferrule size for these canes is 19mm (¾").

PR30089/MARB Marbled black - 36"

PR30089/BLACK Black with tan handle - 36"

PR30089/BAMB Bamboo style - 36"

Max user weight  

2 Comfort Grip Adjustable Crutches

Forearm Crutches have an anatomical soft grip designed to support the palm and spread the pressure, for increasing comfort. Leg and forearm sections adjust independently for a custom fit. Height to handle: 660-890mm (26-35") Handle to cuff: 230-310mm (9-12"). Sold in pairs. Finished in Bronze or Blue metallic. Ferrule size for these canes is 21mm.

PR30081-BZ Bronze

PR30082-BL Blue

Max user weight  

3 Progress Open Cuff Crutch

These lightweight elbow crutches are both stylish and practical. The ergonomically designed handles fit comfortably into the hand and are height adjustable from 720-960mm (28-37"). This crutch also has a height adjustable elbow support 220-290mm (8½-11") Weighs 520g.

Colour: Dark Grey. Sold in pairs.

Ferrule size for these canes is 19mm (¾").

PR30086

Max user weight  

4 Progress Closed Cuff Crutch

These lightweight elbow crutches are both stylish and practical. The ergonomically designed handles fit comfortably into the hand and are height adjustable from 720-960mm (28-37"). This crutch also has a height adjustable elbow support 220-290mm (8½-11") that has a fully pivoting arm cuff. Weighs 520g. Maximum user weight 150kg. Colour: Dark Grey. Sold in pairs.

Ferrule size for these canes is 19mm (¾").

PR30087

Max user weight  





1 Walking Stick Seat – Folding

This useful seat folds flat and clips together to form a walking stick, offering both walking support and a place to sit and rest when needed. Manufactured from aluminium, it is lightweight, strong and corrosion free.

Folded stick height 876mm (34½")

Seat height 571mm (22½")

PR30149

Max user weight  



2 Ferrules

16mm (5/8") - Black.

19mm (3/4") - Black.

19mm (3/4") - Black.

16mm (5/8") - Black. Large Quad Cane Tip

27mm (1 1/8") - Black.

27mm (1 1/8") - Grey.

PR30017-BK

PR30018-BK

PR30034-BK

PR30037-BK

PR30039-BK

PR30039-GR



3 Cane/Crutch Holder

Keeps canes and crutches handy.

A clip-on holder balances the cane and a rubber pad helps keep the holder from sliding off the table. The Cane/Crutch Holder securely grips a cane or crutch with a 16 to 25mm (5/8" to 1") diameter.


PR30158



4 Cane Strap

Cane strap for wood or plastic handled canes

PR30150



1 Deluxe Folding Walker

Folding and lightweight but robust with aluminium frame and adjustable height legs. Folds by use of a simple press down button and becomes only 100mm (4") thick for easy storage or transportation. Silver framed walkers come with plastic handles for easy clean, coloured framed walkers have soft foam, comfort handles. Available in two colours.

Height Adjustment:

Adult Walker 787-965mm (31-38")

Adult Walker 2.4kgs (5.2lbs)

PR30180-BL Blue

PR30180-SL Silver

PR30174/1-SL 5" Fixed Wheels (Pair) - Silver

PR30179/1-BK 5" Swivel Wheels 360° (Pair) - Black

Max user weight  



2 Rigid Walker Frame

Standard rigid walking frame, lightweight but robust with a moulded plastic hand grip and grey high friction rubber tips. Height adjustable from 812-914mm (32-36") with sturdy 25mm (1") anodised aluminium frame. Moulded plastic hand grip. Also available with wheels, fixed or swivel. Width front 457mm (18") Width rear 660mm (26") Depth 546mm (21.5"). Requires assembly.

PR30195

PR30176/A 3" Side Fixed Wheels (Pair)

PR30177/1SL 3" Wheel Extension Legs

PR30173/1 3" Glide Wheels

Max user weight  



3 Hi Riser Folding Frame

Lightweight aluminium walker, with black crossbrace. Can also be used as a toilet safety frame. Width 482mm (19"), height adjusts from 787-914mm (31-36"). Available in silver.

PR30196-SL

PR30174/1-SL 5" Fixed Wheels (Pair) - Silver

PR30179/1-BK 5" Swivel Wheels 360° (Pair) - Black

Max user weight  



1 5" Wheeled Extension Legs

Add 5 inches on to your walker with these fixed, wheeled extension legs, ideal for those who want to add a little height. Fits onto 25mm (1") tube walkers, colour silver, sold in pairs.

Will fit walkers PR30180 - BL and -SL and PR30196- SL

PR30174/1-SL pair

2 360 degree Wheels

360 degree 5" swivel walker wheels, suitable for frames with 25mm (1") tubes. Sold in pairs, colour black.

Will fit walkers PR30180 - BL and -SL and PR30196- SL

PR30179/1-BK pair



3 3" Side fixed wheels for walkers

3" wheeled extension legs. Side-fixed design, suitable for use on the PR30195 Rigid walker, with the tube measuring 22mm ($\frac{7}{8}$ "). Sold in pairs.

Will fit walkers PR30195

PR30176/A pair

4 3" Glide Wheels

These 3" Glide Wheels can make a difference to a walker. Fit on to 25mm (1") tube walkers without any fuss. Colour silver, sold in pairs.

Will fit walkers PR30195

PR30173/1 pair

5 3" Wheel extension legs

3" wheel extension legs to fit on to 25mm (1") tubes. Only available in silver, sold in pairs.

Will fit walkers PR30195

PR30177/1-SL pair



6 Walker Glides

Easy to attach to walkers, fits on rear legs and replaces rubber ferrules. Glides walker across floor. Made from durable plastic.

PR30218-WH 1" White pair

PR30219-GR 1" Grey pair



7 Shopping Trolley

No more struggling up kerbs or steps with heavy old-fashioned shoppers. This new six-wheel shopping trolley has a clever articulated design that just glides up any step or staircase. It's lightweight design (just 2.1Kg) has external pockets for magazines, broly etc., and strap for handbag. Folds flat for easy storage.

PR30281

1 EasyCaddy

The EasyCaddy consists of a tray and an outer case (bucket), which allows you to safely carry personal items and meals easily and with confidence. It is simple to attach to many folding wheeled walkers. It is light sturdy, easy to clean, and it will not upset the balance of the walker.

PR30351



2 The Buckingham Walking Frame Caddy

Developed by healthcare professionals the Buckingham Walking Frame Caddy enables users of wheeled walking frames to safely transport a variety of daily living items. In most instances it can replace the need for a trolley. The caddy enables a patient to carry most items safely and more independently, without being reliant upon carers. The Caddy has one large and one small compartment to allow the user to separate items such as glasses, keys, medication etc. Supplied with a tray, that is designed to safely transport a plate of food, and a mug holder, which catches any spillages. The caddy is easy to fit to most styles of walking frame. Manufactured from a tough, durable plastic it can be cleaned in a dishwasher and has drainage holes, for easy drying. Maximum carrying weight 2kg. Width without tray 434mm (17"). Width with tray 457mm (18"). Depth with tray 310mm (12.25"). Height 200mm (8") Total weight 1kg.

PR30350



3 Synthetic Sheepskin Walker Pads

Cushion hard walker handles.

These soft pads attach to handles easily with a simple hook and loop closure. Machine wash and dry. Sold in pairs.

PR30355



4 Walker Grip Pads

Neoprene foam pads for added comfort on your walker grips. Comes in two different sizes to fit a variety of walkers. Easily removed for washing. Attaches with Velcro.

PR30190 Small Grip

PR30191 Large Grip



5 Walking Frame Net Bag

A knotted nylon bag, that can clip onto any walking frame for carrying shopping or personal belongings.

PR30292





1 Cane & Crutch Holder

Cane & Crutch Holder complete with fixings.

PR30283/CC

2 Bag for Rollator

Optional waterproof, spacious vinyl bag for the Rollator, attaches easily to the frame by straps, Velcro fastening for easy access.

PR30269

3 Basket and Tray for Tri-Walker

Vinyl coated steel basket, which attaches between the handles of the Tri-Wheel. The tray fits neatly on top. Suitable for all Promedics Tri-wheels

PR30353 **Basket**

PR30353/Tray **Tray**



4 Bag For Tri Wheel

A strong and spacious, vinyl waterproof bag optional for the tri wheel walker. Straps around the frame easily, zipper fastening ideal for carrying shopping etc. Suitable for all Promedics Tri wheels.

PR30291




5 Dycem Non-Slip Floor Mat

This non-slip floor mat is non-stick but grips dry, slippery surfaces on both sides. It both grips to the floor and at the same time, prevents the users feet from slipping when standing or when being transferred. Can be used generally in the home in areas where the floor surface is potentially unsafe.

Dycem is non-toxic and can be cleaned with soapy water to retain its properties. It is not effective when wet and must dry before use. All Dycem can be autoclaved at 220°F (103°C).

PR40001 450mm x 600mm (18 x 24") BLUE



Handy Bar

- Security
- Independence
- Mobility
- Confidence
- Injury Prevention
- Emergency Escape
- Universal Fit
- Convenient Size and Weight
- Non Slip Handle
- Patented, Engineered Design

**EXCLUSIVE
TO PROMEDICS**

The customer can get in and out of their vehicle easily and comfortably with this strong simple device. The Handy Bar, with its soft grip, non-slip handle for extra comfort and super strong forged steel shaft, creates a safe, solid support handle where it is needed. The Handy Bar fits the majority of vehicles with the "U" shaped striker. The Handy Bar also incorporates a Seat Belt Cutter and Glass Breaker to help leave a vehicle quickly and easily in the event of an accident

The Handybar® inserts into the U-shaped striker plate on your vehicle door frame. The U-shaped striker plate is one of the strongest parts of the vehicle and is designed to keep your vehicle door closed in the event of a side-impact collision. There is a striker on each of the vehicle's door pillars allowing the Handybar® to be used on both driver and passenger side-doors.

PR40013



Max user weight  





1 Portable Ramps

These wheelchair ramps are designed to be portable, made from anodized aluminium so they are light weight, complete with a carry bag ideal for transporting.

	PR30298/B	PR30299
Length	1.5m (5ft)	2.1m (7ft)
Width	230mm (9")	266mm (10½")
Depth	80mm (3")	50mm (2")
Carrying Wight	6.5kg (14lbs)	15.8kg
Max Load	200kg (29 stone)	400kg (63 stone)
Suitable For:	Wheelchairs Mini Scooters	Wheelchairs All Scooters

Max weight Capacity  




2 Roll Up Wheel Chair Ramp

Aluminium in construction, portable and simple to use quickly makes stairs accessible.

Length	910mm (36½")
Width	750mm (30")

PR30297

Max user weight  



Aids to daily living

 PROMEDICS