



# Orthopaedic

---

Ankle/Back	152
------------	-----

---

Elbow/Hand	153
------------	-----

---

Hand/Wrist	154
------------	-----

---

Knee/Neck/Shoulder	155
--------------------	-----

---

Replay Supports	156-159
-----------------	---------

---

## 1 Anklet

LATEX FREE

A pull on neoprene support for light compression and ankle support. Ideal for ankle sprains, Arthritis, sports injuries and a fractured ankle bone.

		Shoe Size (British)	Shoe Size (Continental)
PR/ANK/02	X Small	3-4	36-37
PR/ANK/03	Small	5-6	38-39
PR/ANK/05	Medium	7-8	41-42
PR/ANK/07	Large	9-10	43-44
PR/ANK/08	X Large	11-12	46-47



## 2 Barnsley Ankle Brace

CONTAINS LATEX

The stitched forefoot section allows this brace to slip on like a sock and the Velcro closure applies variable tension for comfort. Made in quality elastic, giving light ankle support. This brace is especially ideal for ankle sprains, Arthritis, sports injuries and a fractured anklebone.

Code/Left	Code/Right	Size	Shoe Size (British)	Shoe Size (Continental)
PR/EWB/13	PR/EWB/23	Small	3-5	36-38
PR/EWB/15	PR/EWB/25	Medium	5-7	38-41
PR/EWB/17	PR/EWB/27	Large	7-9	41-43



## 3 Pro Tem Belt

CONTAINS LATEX

Simple but effective, the high quality elastic and two flexible posterior stays make this our most popular belt. The edges have quality Lycra binding and a back depth of 20cm (8"), cut away to 18cm (7") at the front for increased comfort when sitting. Especially helps with muscular injuries to the lower back, Lumbago, Osteoporosis, Sciatica and sports injuries.

		Waist Measurements
PR/TB/03	Small	66-74cm (26"-29")
PR/TB/05	Medium	76-84cm (30"-33")
PR/TB/07	Large	86-94cm (34"-37")
PR/TB/08	X Large	96-104cm (38"-41")
PR/TB/09	XX Large	107-114cm (42"-45")



## 4 Sacro Lumbar Support

CONTAINS LATEX

To accommodate the waist/hip differential, two layers of elastic are sewn at an angle giving the belt increased strength and good conformity. The pocket is lined with neoprene for extra comfort and warmth. This belt especially helps with muscular injuries to the lower back, Lumbago, Osteoporosis, Sciatica and sports injuries.

		Waist Measurement
PR/SLS/03	Small	64-74cm (25"-29")
PR/SLS/05	Medium	76-84cm (30"-33")
PR/SLS/07	Large	86-94cm (34"-37")
PR/SLS/08	X Large	96-104cm (38"-41")





**1 Elbow Support**

**LATEX FREE**

This elbow sleeve gives a light support to the elbow area, made from neoprene this sleeve is easy to apply by just pulling on. Especially helps with swollen or painful elbow joints and tennis elbow. Available in three sizes.

		Elbow Circumference
PR/OES/03	Small	18-20cm (7-8")
PR/OES/05	Medium	23-25cm (9-10")
PR/OES/07	Large	28-30cm (11-12")

**2 Tennis Elbow Support**

**LATEX FREE**

This Tennis Elbow support has a high level of patient compliance, as the stretch neoprene accommodates the changing muscle bulk of the moving forearm. With a D-ring strap and adjustable pressure pad, which is easily positioned on the arm to ease the pain of Tennis Elbow. Available in three sizes.

		Circumference 2.5 cm (1") Below Elbow Crease
PR/VPS/03	Small	23-25cm (9-10")
PR/VPS/05	Medium	28-30cm (11-12")
PR/VPS/07	Large	33-35cm (13-14")

**3 Deltaform PLUS Wrist Brace**

**CONTAINS LATEX**

A woven elastic/cotton support, with an aluminium panel inside a vinyl pocket, featuring a longitude seam for excellent conformity to the wrist. The longer Deltaform Plus being 23cm (9") extends to 2/3 of the forearm for radial/ulnac stability. This brace is ideal for broken wrists, forearm fractures, Carpal Tunnel Syndrome, Arthritis, RSI and Tendonitis.

Left	Right	Size In cm	Wrist Circ. in inches	Wrist Circ.
PR/DFP/12	PR/DFP/22	X Small	11.5-13	4 1/2"-5"
PR/DFP/13	PR/DFP/23	Small	14-15	5 1/2"-6"
PR/DFP/15	PR/DFP/25	Medium	16.5-18	6 1/2"-7"
PR/DFP/17	PR/DFP/27	Large	20-21.5	8"-8 1/2"
PR/DFP/18	PR/DFP/28	X Large	23-24	9"-9 1/2"

**4 Wristlet**

**LATEX FREE**

A wrist support made in neoprene, which retains heat and conforms well to give a light compression to soothe painful joints. Wristlet can be cut without fraying to fit swollen joints. Wristlet is ideal for Tendon Injuries, Post Fracture, Arthritis, Scaphoid injury (Wrist injury at base of thumb), Skier's Thumb, Thumb Sprain and RSI.

Left	Right	Size	Circ 2.5cm (1") below elbow crease In cm	Circ. 2.5cm (1") below elbow crease In inches
PR/WL/14	PR/WL/24	Small/Medium	14-18	5 1/2"-6 1/2"
PR/WL/16	PR/WL/26	Medium/Large	17-22	6 1/2"-8 1/2"
PR/WL/LXL/L	PR/WL/LXL/R	Large/X-Large	23-24	8-9 1/2"

**5 Freedom Wrist Brace - Navy**

**CONTAINS LATEX**

A comfortable wrist brace, with a cut away design which allows flexion at the palmar crease and opposition of the thumb to the fingers. Soft Lycra binding prevents chafing of the thumb web. This brace is ideal for broken wrists, forearm fractures, Carpal Tunnel Syndrome, Arthritis, RSI and Tendonitis.

Left	Right	Size	Circ. In cm	Circ. In Inches
PR/FWBN/12	PR/FWBN/22	X Small	11.5-13	4 1/2"-5"
PR/FWBN/13	PR/FWBN/23	Small	14-15	5 1/2"-6"
PR/FWBN/15	PR/FWBN/25	Medium	16.5-18	6 1/2"-7"
PR/FWBN/17	PR/FWBN/27	Large	20-21.5	8-8 1/2"
PR/FWBN/18	PR/FWBN/28	X Large	23-24	9-9 1/2"

## 1 Comfort Cool Thumb CMC Restriction Splint LATEX FREE

Unique strap surrounds the thumb for extra CMC support.

- Made of thin (1.6mm), latex-free, perforated neoprene with a terry cloth liner that helps keep skin cool and comfortable.
- Ideal for CMC joint arthritis and/or instability, tendonitis and thumb discomfort associated with repetitive motion. May be helpful for tenosynovitis of the thumb extensors/abductors.
- Light compression helps reduce thumb pain.
- Contoured strap surrounds the thumb CMC joint and first metacarpal for direct joint support. Allows full hand function to help encourage compliance.
- The thumb sleeve that supports the MP joint can be cut away as needed.
- Hand wash and air dry.
- To size, measure circumference at MP joints.

Left	Right		MP Circumference
PR79562	PR79563	Small	6½" to 7" (17 to 18 cm)
PR79564	PR79565	Med.	7" to 8" (18 to 20 cm)
PR79566	PR79567	Large	8" to 9½"(20 to 24 cm)

## 2 Comfort Cool D-ring Wrist Splint LATEX FREE

Cool, lightweight support provides cumulative trauma relief.

Provides semi flexible support without excessive bulk or warmth. Ideal for wrist disorders such as carpal tunnel syndrome, arthritis and tendonitis. Volar aluminium support is pre-bent to position the wrist in 15° extension and can be removed for adjustments. Circumferentially designed with double layer cotton stockinette under D-ring closures for comfort and easy application. Made of ⅜" (1.6mm) thick perforated neoprene and lined with terry cloth that feels good against the skin. Hand wash and air dry. To size, measure wrist circumference.

Left	Right		Wrist Circumference
PR52962	PR52963	Small	5¾" to 6½" (15-17 cm)
PR52964	PR52965	Medium	6¾" to 7½" (17-19 cm)
PR52966	PR52967	Large	7¾" to 8½" (20-22 cm)

## 3 Comfort Cool Wide Wrist Wrap LATEX FREE

Full coverage for the wrist and distal forearm.

The ⅜" thick perforated neoprene provides support and compression to reduce pain and discomfort. Ideal for soft tissue injuries due to repetitive motion and arthritis. Wide design provides generous overlapping around the wrist and distal forearm. Limits the final degrees of wrist motion while allowing hand function. Neoprene can be trimmed. Soft terry cloth liner feels good against the skin. Hand wash and air dry. Universal sizing.

**PR79510**





1

**1 Standard Knee Sleeve**

**LATEX FREE**

A neoprene sleeve with a reinforced open patella section providing warmth and light compressional support. This sleeve is especially good for Arthritis, mild sprains, after exploratory knee surgery and mild unstable knee.

		Knee Circ. Below Knee	Knee Circumference	Knee Circ. Above Knee
PR/SNO/02	X Small	30cm (12")	30-33cm 12"-13"	38cm (15")
PR/SNO/03	Small	33cm (13")	33-36cm 13"-14"	46cm (18")
PR/SNO/05	Medium	38cm (15")	36-38cm 14"-15"	50cm (19½")
PR/SNO/07	Large	41cm (16")	38-43cm 15"-17"	53cm (21")
PR/SNO/08	X Large	43cm (17")	43-48cm 17"-19"	58cm (23")
PR/SNO/09	XX Large	52cm (20½")	48-56cm 19"-22"	61cm (24")



2

**2 OA Knee Sleeve**

**CONTAINS LATEX**

An elasticated pull on knee sleeve with a soft, cotton - lined foam knee pad, to provide light support and warmth for the arthritic knee. Darts prevent the knee sleeve bunching at the back of the knee. Knee Sleeve used for unstable kneecaps and movement in the kneecaps.

		Knee Circumference
PR/OAS/03	Small	25-28cm (10"-11")
PR/OAS/05	Medium	28-31cm (11"-12")
PR/OAS/07	Large	33-36cm (13"-14")
PR/OAS/08	X Large	38-41cm (15"-16")
PR/OAS/09	XX Large	41-43cm (16"-17")



3

**3 Neck Collar**

**LATEX FREE**

This collar is ideal to stabilise the neck after trauma, or for long term use as a night resting collar. Made in fully washable medium density closed cell foam, which will not collapse with use. Available in four sizes

		Depth x Length
PR/BC/1	Small	7.5 x 43cm (3 x 17")
PR/BC/2	Medium	9 x 47cm (3½" x 18½")
PR/BC/3	Large	9 x 51cm (3½" x 20")
PR/BC/4	Large/Deep	10 x 51cm (4 x 20")



4

**4 Shoulder/Arm Sling**

**LATEX FREE**

This shoulder sling can be used for broken wrist, upper arm and forearm, with its universal design the sling can be adjusted to fit all. A one-piece sling with Velcro closure and stabilising buckle.

PR/UPS      Standard      5cm x 112cm (2" x 44")

## Replay

All Replay supports are made in 'Ventoprene' towelling lined, ventilated neoprene.

The truly ventilated neoprene gives compression and heat retention and allows perspiration to escape, while the absorbent superplush towelling remains comfortable during sporting activity, each support is finished with a stretch Lycra binding for comfort.

The warmth of 'Ventoprene' relieves stiffness and pain, protects joints and muscles, reduces the risk of muscle injury and will retain body heat during intervals of inactivity.

### 1 Replay Lumbar Support

CONTAINS **LATEX**

A lumbar support in Ventoprene providing all round compression and warmth to the sacro-lumbar area. Contoured to the waist for improved fit. This support can be used for muscle injuries to the lower back.

	Size	Waist Circumference in cm	Waist Circumference in inches
PR/RSB/03	Small	61-68 cm	24"-27"
PR/RSB/05	Medium	71-79 cm	28"-31"
PR/RSB/07	Large	81-89 cm	32"-35"
PR/RSB/08	Extra Large	91-99 cm	36"-39"

### 2 Replay Lumbar Support with Panel

CONTAINS **LATEX**

A lumbar support in Ventoprene incorporating a flexible lumbar panel. The two sides pull to provide extra support and compression to the sacro-lumbar area, supporting the spine during flexion, extension and lateral flexion. This support can be used for muscle injuries to the lower back and sacro-lumbar pain.

	Size	Waist Circumference in cm	Waist Circumference in inches
PR/RLB/03	Small	61-68 cm	24"-27"
PR/RLB/05	Medium	71-79 cm	28"-31"
PR/RLB/07	Large	81-89 cm	32"-35"
PR/RLB/08	Extra Large	91-99 cm	36"-39"

### 3 Replay Ankle Support

CONTAINS **LATEX**

This ankle support is made from Ventoprene providing compression and heat retention, incorporating an opening for the heel and contoured at the ankle for improved fit. This ankle support is used for medial/lateral ligament injuries, achilles tendon injuries, malleolus injuries and tendonitis

	Size	Shoe Size
PR/RAS/03	Small	3-6
PR/RAS/05	Medium	6-9
PR/RAS/07	Large	9-12
PR/RAS/08	Extra Large	11-13

### 4 Replay Wrist Support

CONTAINS **LATEX**

An 18cm (7") wrist support in Ventoprene providing compression and heat retention. Contoured at the wrist for improved fit. This support can be used for sprained wrists, twisted wrists, soft tissue injuries, tendonitis wrist protection.

Left	Right	Size	Wrist Circumference in cm	Wrist Circumference In Inches
PR/RWR/13	PR/RWR/23	Small	13-16cm	5"-6½"
PR/RWR/15	PR/RWR/25	Medium	16-21cm	6½"-8½"
PR/RWR/17	PR/RWR/27	Large	20-24 cm	8"-9½"





1

### 1 Replay Calf Support

CONTAINS **LATEX**

This calf support is made from Ventoprene providing compression and heat retention. Contoured at the calf for improved fit. This calf support is used for tennis leg and calf muscle injuries.

	Size	Calf Circumference In cm	Calf Circumference in Inches
PR/RCA/03	Small	30-34 cm	12"-13½"
PR/RCA/05	Medium	34-38 cm	13½"-15"
PR/RCA/07	Large	38-42 cm	15"-16½"
PR/RCA/08	Extra Large	42-46 cm	16½"-18"



2

### 2 Replay Shorts

CONTAINS **LATEX**

Shorts in Ventoprene provide support and warmth during exercise. These shorts can be used for protection of the groin, pelvic and gluteal muscles.

	Size	Waist Circumference in cm	Waist Circumference in Inches
PR/RSS/03	Small	58-66 cm	23"-26"
PR/RSS/05	Medium	68-86 cm	27"-34"
PR/RSS/07	Large	89-102 cm	35"-40"
PR/RSS/08	Extra Large	102-114 cm	40"-45"



3

### 3 Replay Thigh Support

CONTAINS **LATEX**

A thigh support in Ventoprene providing compression and warmth to the quadriceps and hamstring muscles. This support can be used for quadriceps injuries and hamstring injuries.

	Size	Mid Thigh Circumference in cm	Mid Thigh Circumference in inches
PR/RTS/03	Small	33-35 cm	13"-14"
PR/RTS/05	Medium	38-43 cm	15"-17"
PR/RTS/07	Large	46-51 cm	18"-20"
PR/RTS/08	Extra Large	53-58 cm	21"-23"



4

### 4 Replay Long Shorts

CONTAINS **LATEX**

Long shorts in ventoprene provide support and warmth during exercise. These shorts can be used for protection of the groin, pelvic and gluteal muscles and hamstring/quadriceps muscles.


	Size	Waist Circumference in cm	Waist Circumference in inches
PR/RSL/03	Small	58-66 cm	23"-26"
PR/RSL/05	Medium	68-86 cm	27"-34"
PR/RSL/07	Large	89-102 cm	35"-40"
PR/RSL/08	Extra Large	102-114 cm	40"-45"

## 1 Replay Elbow Support

CONTAINS **LATEX**

This elbow support is made from Ventoprene providing compression and heat retention. Contoured at the elbow for improved fit. This elbow support can be used for soft tissue injuries, throwers elbow, sprained elbow and elbow protection.

	Size	Forearm Circumference in cm	Forearm Circumference in inches
PR/RES/03	Small	23-25 cm	9"-10"
PR/RES/05	Medium	25-30 cm	10"-12"
PR/RES/07	Large	30-35 cm	12"-14"
PR/RES/08	Extra Large	35-41 cm	14"-16"



## 2 Replay Padded Elbow Support

CONTAINS **LATEX**

A padded elbow support in Ventoprene providing compression and heat retention. Padding at the elbow reduces the effect of friction burns on court, fields or Astroturf and protects against knocks, bruising and abrasion. This support can be used for throwers elbow, sprained elbow, elbow protection and soft tissue injury.

	Size	Forearm Circumference in cm	Forearm Circumference in inches
PR/RPE/03	Small	23-25 cm	9"-10"
PR/RPE/05	Medium	25-30 cm	10"-12"
PR/RPE/07	Large	30-35 cm	12"-14"
PR/RPE/08	Extra Large	35-41 cm	14"-16"



## 3 Replay Tennis Elbow Support

CONTAINS **LATEX**

A tennis elbow support in Ventoprene providing compression and heat retention. The adjustable pressure pad can be positioned to ease the pain of medial or lateral epicondylitis. Stretch Ventoprene accommodates the changing muscle bulk of the moving forearm when playing racquet sports. This support can be used for tennis elbow, golfers elbow and elbow protection.

	Size	Forearm Circumference in cm	Forearm Circumference in inches
PR/RVP/03	Small	23-25 cm	9"-10"
PR/RVP/05	Medium	25-30 cm	10"-12"
PR/RVP/07	Large	30-35 cm	12"-14"
PR/RVP/08	Extra Large	35-41 cm	14"-16"



## 4 Replay Hinged Knee Support

CONTAINS **LATEX**

A hinged knee support in Ventoprene, with metal hinges to control medial and lateral instability, and a patella opening to stabilise tracking. Velcro straps above and below the knee provide improved fit. This support can be used for unstable medial/lateral ligaments, cartilage damage, runners knee and knee protection.

	Size	Mid Patella Circumference in cm	Mid Patella Circumference in inches
PR/RHN/03	Small	30-34 cm	12"-13½"
PR/RHN/05	Medium	34-38 cm	13½"-15"
PR/RHN/07	Large	38-42 cm	15"-16½"
PR/RHN/08	Extra Large	42-46 cm	16½"-18"







1

**1 Replay Padded Knee Support** CONTAINS **LATEX**

A padded knee support in Ventoprene providing compression and heat retention. Padding at the knee reduces the effect of friction burns on court, field or Astroturf and protects against knocks, bruising and abrasions. This support can be used for knee protection, cartilage damage, soft tissue injuries and tendonitis.

	Size	Mid Patella Circumference in cm	Mid Patella Circumference in inches
PR/RNC/03	Small	30-35 cm	12"-14"
PR/RNC/05	Medium	35-38 cm	14"-15"
PR/RNC/07	Large	38-43 cm	15"-17"
PR/RNC/08	Extra Large	43-48 cm	17"-19"



2

**2 Replay Open Knee Support** CONTAINS **LATEX**

An open knee support in Ventoprene providing compression and heat retention. Use for soft tissue injury, runners knee.

	Size	Mid Patella Circumference in cm	Mid Patella Circumference in inches
PR/RNO/03	Small	30-35 cm	12"-14"
PR/RNO/05	Medium	35-38 cm	14"-15"
PR/RNO/07	Large	38-43 cm	15"-17"
PR/RNO/08	Extra Large	43-48 cm	17"-19"



3

**3 Replay Stabilised Knee Support** CONTAINS **LATEX**

A stabilised knee support in Ventoprene with spiral stays to give medial and lateral support, and a patella opening to stabilise tracking. This support can be used for unstable medial/lateral ligaments, cartilage damage, runners knee and knee protection.

	Size	Mid Patella Circumference in cm	Mid Patella Circumference in inches
PR/RST/03	Small	30-35 cm	12"-14"
PR/RST/05	Medium	35-38 cm	14"-15"
PR/RST/07	Large	38-43 cm	15"-17"
PR/RST/08	Extra Large	43-48 cm	17"-19"




4

**4 Replay Knee Support** CONTAINS **LATEX**

A knee support in Ventoprene providing compression and heat retention, contoured at the knee for improved fit. This support can be used for cartilage damage, soft tissue injuries, knee protection and tendonitis.

	Size	Mid Patella Circumference in cm	Mid Patella Circumference in inches
PR/RSN/03	Small	30-35 cm	12"-14"
PR/RSN/05	Medium	35-38 cm	14"-15"
PR/RSN/07	Large	38-43 cm	15"-17"
PR/RSN/08	Extra Large	43-48 cm	17"-19"



Aids to daily living

 PROMEDICS