



Eating & Drinking Assistance

Cutlery	92-94
---------	-------

Bibs & Eating	95
---------------	----

Plates	96-97
--------	-------

Cups & Beakers	97-102
----------------	--------

1 Good Grips Cutlery

Good Grip cutlery has a soft cushion grip that keeps the utensil in the hand - even when wet!

The built up handles are made of soft latex - free material, with flexible ribbing that's comfortable and adapts to any grip.

Each stainless steel utensil has a special twist built into the metal shaft that allows the spoons and forks to bend to any angle. You can bend them for either Left or Right handed use. The rocker knife requires only minimal arm strength for cutting. A special extended safety cap protects fingers. Dishwasher Safe.

- PR65591 Fork each**
- PR65595 Youth Spoon each**
- PR65593 Teaspoon each**
- PR65594 Desert Spoon each**
- PR65592 Rocker Knife each**
- PR65592/S Serated RockerKnife each**



Good Grips Weighted Cutlery

Textured handle is secure in the hand. Each cutlery piece features the Good Grips handle and a special twist built into the metal shaft allows it to be bent for either right or left handed use. The angled rocker knife requires only minimal arm strength for cutting. A special extended safety cap protects fingers. The 6oz. (170g) of weight in the built-up-handle provides more control so that the utensil reaches the mouth easier. Ideal for persons with limited hand control, Parkinsons disease or spasticity.

- PR65560 Weighted Fork**
- PR65561 Weighted Teaspoon**
- PR65562 Weighted Desert Spoon**
- PR65564 Weighted Souper Spoon**
- PR65565 Weighted Rocker knife**

2 Good Grips Coated Spoons

Idea for persons with biting reflexes and limited co-ordination. Good grips coated spoons have built up handles that are easy to hold even when wet! Plastisol coating protects teeth and lips. Ideal for clients with spasticity or limited hand control. Not recommended for those with heavy biting reflexes. Can be used with the Goodie-Strap (below). Hand washing in cold water prolongs product life.

- PR65596 Youthspoon**
- PR65597 Teaspoon**
- PR65598 Desert Spoon**
- PR65600 Sample Kit - Contains:**
1 x PR65596, 1 x PR65597, 1 x PR65598

3 Goodie-Strap

Fits large and built-up handles. This elastic strap has two different-sized loops that fit various handle sizes. The regular strap fits most adult hands and will accommodate Good Grips utensils and other large, built-up handles. The paediatric strap fits the Good Grips Paediatric Spoon and other utensils with an approximate circumference of 40 to 76mm (2½" to (3")), and a handle length of at least 76mm (3"). Machine wash and air dry.

- PR65599 Regular**
- PR65589 Paediatric**

4 Good Grips Paediatric Spoon

Specially sized for children. This paediatric spoon has the same features as the Good Grips Spoons, with a smaller handle and bowl. Institutional dishwasher safe. Total length measures 6" (150mm). Can be used with the Paediatric Goodie-Strap PR65589

- PR65588 Paediatric Spoon**



5 Good Grips Souper Spoon

Raised lip is ideal for persons with tremors or poor hand control.

The innovative lip of the Souper Spoon is specially designed to prevent spills while moving the spoon from the plate to the mouth. Liquids and food stay on the spoon where they should be! The one-piece, stainless steel lip is easy to clean and fits in the mouth.

The Good Grips Souper Spoon features the Good Grips handle and the special twist in the metal shaft. The total length measures 180mm (7½"). The Good Grips model is ideal for persons with Alzheimers, CVA, Parkinsons disease, arthritis, or neurological impairments.

Înstitutional dishwasher safe.

- PR65590**



1 Comfort Grip Cutlery

Nonslip grips rotate for better control

These utensils are ideal for anyone with a poor, weak grip, allowing maximum control with minimum effort, giving independence at mealtimes. Each utensil is lightweight, 85grams (3oz) each, with a sturdy, non-slip, 7/8" (22mm) diameter grip handle which rotates to accommodate most eating difficulties. The ridges on the handle help to position the fingers for a more secure grip. The blade of the Comfort Knife has been specially designed to be useful, both as a regular dinner knife and as a rocker knife. It has a curved cutting edge making it a suitable cutting knife for a one-handed person. Not only suitable for use in rehabilitation centres, but also attractive enough to use in the home. These utensils are institutional dishwasher safe and Latex free. Each utensil can be purchased separately.

PR65631 Straight Fork

PR65632 Straight Knife

PR65633 Straight Teaspoon

PR65634 Straight Desert Spoon

PR65635 Left-Handed Fork

PR65636 Left-Handed Teaspoon

PR65637 Left-Handed Dessertspoon

PR65638 Right-Handed Dessertspoon

PR65639 Right-Handed Teaspoon

PR65640 Right-Handed Fork

2 Soft Handled Cutlery

These soft handled utensils have large grips made of closed cell foam which will not absorb moisture. Lightweight, and ideal for people with arthritis or tender hands and finger joints. The foam grip slides off easily for cleaning.

PR65641 Soft Handled Knife

PR65642 Soft Handled Fork

PR65643 Soft Handled Teaspoon

3 Lifestyle cutlery

Cutlery that truly fits your hand. Perfect for anyone with arthritis or some limited hand grasp, yet attractive enough for family use.

This soft, comfortable to the hand cutlery has been designed to stay in place when wet or when turning your grip. The ergonomically designed contours allow for a comfortable fit in almost any hand while adding to stability. The slight curve at the end allows the grip to nestle securely in the hand giving a feeling of security. Available in grey. Pieces sold separately.

PR65655/K Knife

PR65655/F Fork

PR65655/SS Soup Spoon

PR65655/TS Tea Spoon

1 Rocking "T" Knife

Easy cutting for people with a weak grasp. This stable, easy-to-use utensil cuts by rocking the knife on food. Comfortable wooden handle is shaped to fit the hand. The safety edge, Stainless steel blade may be sharpened with a Knife sharpener. Household dishwasher safe.

PR65021

2 Plastic Coated Spoons

Coating protects teeth and lips. The special plastisol coating on these plastic coated spoons protects teeth and lips during use. Excellent for clients with Spasticity or limited hand control. Not recommended for those with heavy biting reflexes. Handle measures 13mm ($\frac{1}{2}$ ") in diameter at base. Hand washing in cold water will prolong product life.

PR65628 Youthspoon

PR65629 Teaspoon

3 Care Spoons

Small spoons protect teeth. Durable, lightweight spoons provide exceptional protection for teeth and gums during feeding. Small bowls helps control food intake. Spoons measure 150mm (6") long and 13mm ($\frac{1}{2}$ ") in diameter. The small spoon has a 25mm (1") bowl and the large spoon has a 32mm ($1\frac{1}{4}$ ") bowl. Household dishwasher safe.

PR65010 Small (10)

PR65011 Large (10)

4 Hardcoated Spoons

Helps protect teeth and lips while eating. The plastic coating on these hardcoated spoons protects teeth and lips of clients with biting reflex, limited hand control or hypersensitivity to temperature. Not recommended for person with heavy bite reflexes. Sturdy plastic material keeps edges smooth, even after washing. Handle measures $\frac{1}{2}$ " (13mm) in diameter at the base. Hand washing in cold water will prolong product life.

Spoon	Description	Code
A	Dessert Spoon	PR65615
B	Long Spoon	PR65613
C	Tea Spoon	PR65612
D	Soup Spoon	PR65614

5 SteadySpoon

Spoon bowl swivels to keep food stable and level despite tremors or weakness.

An active counterbalance and weight steadies the spoon bowl, keeping food level and in place. Using the SteadySpoon helps promote productive, independent feeding. Includes a soft Velcro strap that secures the SteadySpoon to the hand and can be adjusted as needed. Ideal for anyone with a weak grasp, hand tremors or limited range of motion. The SteadySpoon weighs 227g (8oz). The handle measures 150mm long and 32mm in diameter ($5\frac{3}{4}$ " x $1\frac{1}{4}$ "), and the length including spoon is 250mm (10"). Spoon is top rack dishwasher safe.

Paediatric version for age 4 and under, is slightly smaller, handle measures 111mm ($4\frac{1}{2}$ ") total length including spoon 200mm ($7\frac{3}{4}$ ") and weighs 5oz.

PR65000 SteadySpoon adult


PR25020 SteadySpoon paediatric


6 Unbreakable Spoons

Tough Nylon material is kind to lips and teeth- durable enough to be used by individuals with strong bite reflexes. Utensils can be customised for user's need by cutting to shape with scissors and smoothing with an emery board. Can be cleaned in a dishwasher up to 180°F and sterilised. Scissors and emery board are not included. Latex Free. Package of 10.

PR65647 Small Plastic Spoon for bite reflex PK 10

PR65648 Large Plastic Spoon for bite reflex PK 10





1

1 Universal Cuff

Economical, adjustable holder
Cutlery, pens and pencils can be held with the Utensil Holder. Hook & loop closure with a D-ring strap adjustment for palms measuring 180-200mm (7 - 8") in circumference. Easy to put on and take off. Pocket measures 76mm (3") long and is sewn closed on one end. Made of comfortable cotton. Machine wash, gentle cycle.

PR65058



2

2 Coloured Foam Tubing

Bright colour assortment is ideal for visually impaired and paediatric use.

Use Coloured Foam Tubing to build up cutlery and other hard-to-grasp handles. This closed-cell, slip resistant foam is available in three different diameters shown in the chart below. All tubing is sold in a package of six 300mm (12") long tubes. Household, dishwasher safe.

	Colour		Outside Diameter	Hole Size
PR65007	Tan	(6)	1" (2.5cm)	¼" (6.4mm)
PR65008	Red	(6)	1 ⅛" (2.9cm)	⅜" (9.5mm)
PR65009	Blue	(6)	1 ⅛" (2.9cm)	¾" (19mm)
PR65069	Assorted	(6)	Assorted	Assorted



3

3 Dinner Bibs

These colourful bibs will brighten mealtimes with a variety of patterns and assorted colours. They measure 460mm (18") wide x 890mm (35") long, with a poly/cotton front and a waterproof backing. The neck of these bibs open and close with snap buttons. Colours may vary.

PR65645



4

4 Terry Cloth Bib

Large, soft bib allows for neater, more comfortable mealtimes. Hook and loop closures makes this bib easy to put on. Made of 90% cotton 10% polyester. Reinforced with cotton binding. Bib measures 450 x 760mm (18" x 30"). Machine washable. Meets compliance law standards.

PR65001

5 Terry Cloth/Vinyl Bib

Provides protection from spilled liquids. The reversible Terry Cloth/Vinyl Bib has terry cloth on one side and soft vinyl on the other. Neck opening fasteners with hook and loop closures. Bib measures 18" x 27" (450 x 680mm). Machine washable. Meets compliance law standards.

PR65002



5

1 GripWare Round Scoop Dish

The low front and high back are ideal for persons with limited motor co-ordination or the use of only one hand.

Dish measures 200mm (8") in diameter. Available in plastic or melamine. The plastic dish is microwave safe and can be gas autoclaved. The melamine dish should not be microwaved, and is safe up to 125°F (52°C) in the dishwasher.

PR65006 Plastic

PR65016 Melamine

2 Scoop Plate

Rim with grooved lip to help hold the utensil.

Dish measures 230mm (9") in diameter. The base has non-slip rubber feet that hold the dish in place. Made of melamine, do not microwave. Dishwasher safe up to 125°F (52°C).

PR65627

3 Partitioned Scoop Plate

Compartments help to keep food separate and provide more surfaces for scooping.

Dish measures 220mm (8½") in diameter with 19mm (¾") high dividers. Two 118ml (4fl.oz.) compartments and one 237ml (8fl.oz.) compartment. Institutional dishwasher safe. The Scoop dish is made of melamine, do not microwave. Dishwasher safe up to 125°F (52°C).

PR65601

4 Inner Lip Plate and Suction Lip Plate

Encourage independent eating.

Special lip design prevents spills by letting food be pushed to the side of the plate, then scooped up with a fork or spoon. Ideal for one-handed dining. Both plates have a depth of 13mm (½") and a diameter of 230mm (9"). The Inner Lip Plate is institutional dishwasher safe and may be microwaved or autoclaved. The Suction Lip Plate secures to the table using suction feet, is household dishwasher safe and should not be microwaved.

PR65017 Inner Lip Plate

PR65019 Suction Lip Plate

PR65026 Pack of 3 Replacement Suction Cups

5 Scooper bowl

Rolled edge helps children eat without spills.

Prevent spills and overflow with the scooper bowl. Children will love the unique shape. Rubber - ring base prevents sliding. Bowl has a 130mm (5") diameter. Unbreakable. Top rack household dishwasher safe.

PR65004

6 Clear Food Guard

Attractive enough for any table.

The inconspicuous design for the Clear Food Guard makes it ideal for home use or when dining out in a restaurant. Clear polycarbonate plastic is strong and functional, yet attractive. A perfect tool for persons who need a little extra assistance when eating.

A snug fit means that the Clear Food Guard stays on the plate securely when food is scooped up against it so that no food slips off the plate. Tapered ends prevent the food guard from interfering with utensils. The Food Guard can be attached quickly even by clients with the use of just one hand. The grooved edge simply slides onto the plate, forming a gapless fit. Fits plates ranging from 210 to 250mm (8½" to 10") in diameter. Designed to nest together, taking up less storage space. Top rack, household dishwasher safe. Sold individually or in a package of 10.

PR65029 Individual

PR65574 Pack of 10





1

1 Food Guard

This flexible food guard can be fitted to ordinary plates to assist with one handed eating. Three clips attach it to the rim of the plate and create a high inward sloping face to help with food collection and prevent spillage. Colour White only. Plate size 190 to 254mm (7½ to 10") external diameter. Microwave and dishwasher safe. Weight 45g.

PR65071



2

2 Stay Warm Feeding Dish

Keep food warm for longer using this convenient dish. A chamber under the plate is easily filled with warm water, keeping food warm for longer. The dish also has a non slip ring which keeps the dish in place, reducing the risk of spillage. Do not fill the plate with boiling water.

PR65110-N



3

3 Melamine Dish

This melamine dish measures 8" in diameter and is designed to help with food collection. The dish is not suitable for microwaves, but is dishwasher safe.

PR65653



4

4 Teapot Beaker

This white plastic feeding cup with two handles has a partially covered lid to prevent spills. The teapot style spout enables the cup to be used at a less steep angle when drinking.

Capacity 250ml (9fl.oz).

PR65122



5

5 Two Handled Mug

Each mug has two handles to assist control, and is supplied with two lids, one for anti-spill and one with a drinking spout. The mug is dishwasher safe, but not the lid. Colours may vary. Capacity 250ml (9fl.oz). Weight 110g.

PR65123

6 Clear Mug with Handle

A tough, clear plastic mug, with large, open handle to grip round or under without making a fist. This mug is supplied with two lids, one for anti-spill and one with a drinking spout. It is of a simple design, stackable and easy to clean. The mug is dishwasher safe however the lids are not, but can be used in the microwave. Capacity 300ml.

PR65100



6

1 Plastic Cup With Lid

Choose the lid that suits your clients specific needs. The sturdy, textured Plastic Cup with Lid is available with either a Snorkel Lid or Snap On Lid. Holds 237ml (8fl.oz.). Replacement lids are available in a package of 6. Institutional dishwasher safe.

PR65054 Cup w/ Snorkel Lid

PR65035 Snap On Lids (6)

2 Nose Cut-Out Beaker

Unique design enables drinking without tilting the head. The Nose Cut-Out Beaker is ideal for people who choke easily or cannot tilt their heads back due to neck injuries. Useful for children who have cerebral palsy. The cut-out fits around the nose when drinking. The 73mm (2 $\frac{7}{8}$ "") wide opening allows the therapist to observe the liquid going into the clients mouth. Available in a regular or clear model. The regular model is made of unbreakable plastic. The clear model is made of polypropylene plastic. Holds 237ml (8fl.oz.). Household dishwasher and microwave safe.

Regular

PR65051 Single Beaker

PR65052 Pkg. of 10

PR65067 Clear

3 Flexi Cut Cup

Drink without neck extension.

Unlike other cut cups, the Flexi Cut Cup can be squeezed gently to change the shape of the cup lid. Made of flexible plastic with a cutout design. Top rack household dishwasher safe. Sold in a package of 5.

PR65055 Pink 1 fl. oz. (30ml)

PR65053 Blue 2 fl. oz. (59ml)

4 Spill Proof Kennedy Cup

Designed specifically for rehabilitation use.

This lightweight, easy-to-grip Spill Proof Kennedy Cup has a sturdy, screw-on lid. This cup will not spill and can be used in either a reclined or seated position. Liquid level is visible and can be monitored easily. Holds up to 207ml (7fl.oz.). Household dishwasher safe and autoclavable.

PR65025

5 Feeding Cup With Flow Control


Regulate drink flow and intake.

The soft rubber button release on the Feeding Cup allows control of liquid flow and intake. Cup holds 237ml (8fl.oz.). Replacement buttons are available in a package of 12. Institutional dishwasher safe and autoclavable.

PR65045 Cup

PR65046 Buttons (12)





1 Feeding Cup and Base

Sturdy base prevents tipping.

Prevent dribbles and spills with the Feeding Cup and Base. The mouthpiece hole will accommodate a drinking straw. The cup is securely positioned with a 130mm (5") diameter base. Cup holds 237ml (8fl.oz). Made of sturdy polyethylene plastic. Cup and base are sold separately. Institutional dishwasher safe and autoclavable.

PR65043 Cup

PR65044 Base

2 Dysphagia Mug

Unique shape provides control for persons with dysphagia.

This mug is designed to direct the flow of liquid to the centre of the mouth. Allows for drinking liquids or thickened liquids without tilting the head back. Large handle is easy to hold. Hand easily fits through the handle, making it ideal for persons with weak grasp. Holds 237ml (8fl.oz.). Weighted base adds stability. Household dishwasher and microwave safe. Autoclavable. Colour may vary.

PR65047

3 Doidy Cup

Unique angle reduces the need to tilt the head when drinking.

Children will love drinking from the Doidy Cup. The unique shape is appealing to children while making drinking easy. Two handles fit small hands. Holds up to 192ml (6½fl.oz.). Household dishwasher safe. Colour may vary.

PR65018

4 Drinking straws

Adaptable straws are reusable.

These versatile drinking straws are available in large and small Plexiglas, and flexible polyethylene. All straws measure 460mm (18"). The small straw has a 3.2mm (¼") hole and the large straw has a 6.4mm (¼") hole. The Plexiglas can be adjusted to any angle using a heat gun. The flexible straws have a 4.95mm (3/16") hole and can be shaped easily by hand.

PR65040 Large (5)

PR65039 Flexible (10)

5 One-Way Straws

Reduce air intake with these straws.

Liquid will not flow back down these straws after the user removes lips. The one-way valve eliminates sucking air from an empty straw. Ideal for persons with Parkinson's disease. Attach to glasses and mugs with the included clips. Each package includes one 250mm (10") and one 170mm (7") straw.

PR65059

6 Novo Cup

Perfect for people who want to drink without having to be helped.

The Novo Cup allows you to drink comfortably while lying down.

Your are in complete control of the rate of flow and don't need to be propped up by pillows or helped by carers. Ideal for those with restricted movement of the head and neck.

- Spill-proof when tipped over, so fewer changes of bed linen.
- Takes both hot and cold liquids
- Transparent so contents are visible
- Measurement markings
- Holds 250ml of liquid
- Dishwasher, Microwave and autoclave safe

PR65075

1 Thumbs Up Cup

Ideal for people with arthritic hands

Two thumb rests make this a uniquely designed cup ideal for the arthritic hand. Insulated, moulded plastic keeps liquids hot or cold longer, while keeping the outside of the cup comfortable to the touch. Includes lid. Holds 325ml (11fl.oz). Weighs 118ml (4oz). Household dishwasher and microwave safe.

PR65070

2 Sure Grip Cup

This stackable, polycarbonate mug features deep grooves allowing a comfortable grip in the hand and resists slippage. It is suitable for all age groups, hot and cold liquids and the virtually spill-proof, patented valve in the lid prevents spillage. Only a small drop escapes before the vacuum seals off the main flow. The Sure Grip Mug has a spout, angled to facilitate easy use. Another feature of this cup is the inclusion of a 100ml and 200ml measure guide on the outside, this comes in useful when volumes of liquid need to be recorded. Microwave and dishwasher safe, available in three colours. The optional lid changes colour if the contents are in excess of 50°C. in temperature.

PR65571 Clear

PR65572 Blue

PR65573 Pink

PR65580 Temperature regulated lid

3 Plastic drinking straw holders

Spring clip fits any size cup.

Snap this straw holder to the edge of any glass. Designed to hold all straws of all sizes at comfortable angles. The plastic drinking straw holder has eight holes and attaches easily with a spring clamp fastener.

PR65041 Plastic

4 Two Handled Mug

Monitor liquid levels with this see-through mug.

Use this clear plastic mug with or without the included lid. The spout in the lid directs the flow of liquid. The two large handles provide for an easy grasp and the wide base helps prevent accidental tipping. Holds 10 fl. oz. (295ml). Household dishwasher and microwave safe. Replacement lids sold separately.

PR65616 Mug

PR65617 Replacement Lid

5 Handycup

The rim of the HANDYCUP is angled so that a person can drink without tipping their head back to drain the glass. It can be tipped beyond the point where a normal cup would hit the bridge of your nose. The cup is transparent so that the level of the contents is always visible. The cup has large handles so that the person can fit their hand inside for better grip or can hold handles. The large base helps to prevent tipping. Designed for people with a poor grip, limited neck movement, a bed patient, or a person in a wheelchair with a head-rest.

PR65646 Clear

PR65646/B Blue

PR65646/G Green

PR65646/O Orange

PR65646/LID Clear lid





1

1 Melaware Cup

Persons with limited grasp can hold this cup confidently. The thumb and index finger slip easily around the narrow stem of the Melaware Cup. Ideal when finger strength is decreased. Holds 237ml (8fl.oz.). Melamine cup is institutional dishwasher safe but cannot be microwaved.

PR65024



2

2 Insulated Mug

Keep drinks at the right temperature. The Insulated Mug With Lid helps keep liquids hot or cold longer. Large handle makes grasping easier and the sealing lid helps prevent spilling. Straw hole directs the flow of liquid. Made of durable polypropylene. Holds 142ml (12fl.oz.). Household dishwasher and microwave safe.

The Weighted Insulated Mug With Lid has the same features as the insulated Mug With Lid, with the added feature of a weighted rubber base that helps prevent accidental tipping. Ideal for people with tremors. Weighs 269g (9½oz.). Household dishwasher safe but cannot be microwaved. Replacement lids for both mugs are available in a package of 6.

PR65030 Insulated Mug

PR65630 Weighted Insulated Mug

PR65035 Replacement Snap On Lids (6)

PR65068 Snorkel lids (6)



3

3 Insulated Mug

A quality insulated mug at a great value. This polypropylene mug keeps liquids hot or cold as necessary. Foam insulated and break-resistant. Holds 237ml (8fl.oz.). Can be used with the Snorkel or the Snap On Lid, sold separately in a package of 6. Household dishwasher and microwave safe.

PR65031 Insulated Mug

PR65035 Snap On Lids (6)



4

4 Feeding Cup With Long Spout

Large spout for intake without dribbling. The Feeding Cup with Long Spout features a mouthpiece that is large enough to hold a straw. The spout can be cut down as skills improve. Place a finger over the air vent to easily regulate the flow of liquid. Holds 237ml (8fl.oz.). Made of polypropylene plastic. Available with or without handles. Replacement lids are sold in a package of 3. Household dishwasher safe.

PR65033 Cup

PR65032 Cup with handle and 2 lids narrow and wide spouts

PR65033-3L Lids (3)

PR65034 Cups (20)



1 China Cup and Saucer

This lightweight attractive bone china cup and saucer is suitable for those with poor grip. The cup has handles on both sides to provide additional stability. Microwave and Dishwasher safe.

Weight Cup: 190g (6.7 oz) Saucer: 130g (4½ oz)
Volume: 160ml

PR65220



2 Easi 2 Drink

Difficulty with drinking? Problems with tremors? This innovative device has been specially designed to aid drinking without spilling. A plastic insert placed in the mug restricts wave motion and the tendency for the liquid to spill. Inserts are manufactured in food quality plastic and will withstand 0-100°C temperature range. Inserts are easily removed for cleaning and can be used over and over again. Additives such as sugar and milk can be added and stirred before pushing the insert into the mug. Each ceramic mug comes complete with insert and stainless steel extractor. The mug and insert are dishwasher safe.

PR65221



3 Nosey Cup

The rim of the Nosey Cup is angled so that a person can drink without tipping their head back to drain the glass. It can be tipped beyond the point where a normal cup would hit the bridge of your nose. The glass is transparent so that the level of the contents is always visible. The bottom of the bowl is rounded to assist draining the contents and to make washing easy. The short stem and tapered cup make the glass easy to grasp. The large base helps to prevent tipping. Designed for people with a poor grip, limited neck movement, a bed patient, or a person in a wheelchair with a head-rest. Capacity 8oz (260ml).

PR65644

